

SORT OUT IX

A randomized trial comparing BioFreedom with Orsiro® in an all-comers patient population

Conclusions

- BioFreedom* was not non-inferior to Orsiro in the primary endpoint of Target Lesion Failure (TLF) at 12-month follow-up (5.2% vs. 4.0%, p-non inferiority = 0.123)¹, and no significant difference was found at 24-month follow-up (6.3% vs. 7.8%, RR 1.23 95% CI 0.94 – 1.61)².
- At 2 years, Target Lesion Revascularization (TLR) rate was significantly lower in the Orsiro stent group compared to the BioFreedom stent group (2.6% vs. 5.1%, RR 1.98 95% CI 1.36 – 2.89)².
- BioFreedom and Orsiro both had similar safety and risk profile for definite Stent Thrombosis (ST) up to 2 years.²

Study design

Randomized, multi-center, single blind, all-comers, two-arm, non-inferiority trial comparing BioFreedom to Orsiro stent in patients treated with PCI at 4 hospitals in Denmark.

Endpoints


Primary endpoint

Target Lesion Failure (TLF) defined as the composite of:

- Cardiac Death
- Myocardial Infarction (MI) not related to any segment other than the target lesion
- Target Lesion Revascularization (TLR)

Secondary endpoints

- Individual components of the primary endpoint
- Stent Thrombosis (ST) rate according to the ARC definition^Δ
- All-case death
- Target Vessel Revascularization (TVR)
- POCE (Death, MI, or any revascularization)





Patient characteristics ¹	BioFreedom n = 1,572	Orsiro n = 1,579
Age (years)**	66.4 ± 10.7	66.1 ± 11.1
Male	77.5%	77.3%
Diabetes	19.3%	19.2%
Current smoker	29.8%	29.3%
Prior PCI	20.9%	20.9%
Prior CABG	8.4%	7.0%
Prior MI	14.7%	15.2%
Stable angina	42.7%	40.8%
NSTEMI / Unstable Angina	28.9%	28.7%
STEMI	23.3%	25.1%
Other	5.1%	5.3%

Lesion and Procedural characteristics ¹	BioFreedom n = 1,966 [◇]	Orsiro n = 1,985 [◇]
Lesion per patient**	1.3 ± 0.6	1.3 ± 0.6
Lesion type B2/C	60.6%	58.1%
Reference vessel size (mm)**	3.3 ± 0.6	3.3 ± 0.6
Number of stents		
Per patient**	1.6 ± 0.9	1.6 ± 0.9
Per lesion**	1.3 ± 0.6	1.2 ± 0.6
Total stent length (mm)		
Per patient**	31.1 ± 21.9	30.6 ± 19.8
Per lesion**	24.7 ± 16.0	24.3 ± 13.6

* BioFreedom is a trademark or registered trademark of Biosensors International Group, Ltd.
 ** Data shown as mean – SD
 ◇ Number of lesions
 Δ According to Academic Research Consortium (ARC) criteria for acute, subacute, late, very late and cumulative stent thrombosis




TLF at 12- and 24-month²



TLF Components at 24-month


	BioFreedom n = 1,570	Orsiro n = 1,576	Rate ratio (95% CI)
Cardiac Death	2.0%	2.6%	0.78 [0.49 - 1.24]
MI	4.3%	4.4%	0.97 [0.69 - 1.36]
TLR	5.1%	2.6%	1.98 [1.36 - 2.89]

TLR 24-month²



ST 24-month²

RR 0.80
95% CI 0.47-1.37



BioFreedom



Orsiro

Subgroup²

		Rate of TLF events BioFreedom Orsiro	Rate Ratio (95% BC ¹ **)	p for interaction
ACS	no	63 (8.4%) 48 (6.6%)	1.29 (0.88-1.88)	0.75
	yes	59 (7.2%) 52 (6.1%)	1.18 (0.81-1.71)	
Age	≤ 65	32 (4.8%) 33 (4.6%)	1.06 (0.65-1.73)	0.53
	> 65	90 (9.9%) 67 (7.8%)	1.28 (0.93-1.76)	
Diabetes mellitus	no	83 (6.5%) 67 (5.3%)	1.26 (0.91-1.74)	0.83
	yes	39 (12.8%) 33 (10.9%)	1.17 (0.73-1.88)	
LAD	no	66 (8.4%) 48 (6.1%)	1.39 (0.96-2.02)	0.36
	yes	56 (7.1%) 52 (6.6%)	1.09 (0.74-1.59)	
Lesion type C	yes	60 (10.1%) 50 (8.5%)	1.18 (0.81-1.72)	0.72
	no	62 (6.4%) 49 (5.0%)	1.30 (0.89-1.90)	
Male	no	23 (6.5%) 24 (6.7%)	0.97 (0.54-1.72)	0.35
	yes	99 (8.1%) 76 (6.2%)	1.32 (0.97-1.78)	
MVD	no	95 (7.3%) 77 (5.9%)	1.25 (0.92-1.69)	0.88
	yes	27 (10.3%) 23 (8.6%)	1.19 (0.68-2.09)	
One stent per patient	no	75 (7.5%) 61 (6.1%)	1.25 (0.89-1.75)	0.83
	yes	44 (7.8%) 38 (6.7%)	1.18 (0.76-1.82)	
Previous MI	no	90 (6.9%) 75 (5.8%)	1.21 (0.89-1.65)	0.45
	yes	29 (12.9%) 21 (9.0%)	1.45 (0.83-2.54)	
Previous PCI	no	80 (6.6%) 71 (5.7%)	1.15 (0.83-1.58)	0.66
	yes	39 (12.1%) 25 (8.0%)	1.53 (0.92-2.54)	
STEMI	no	104 (8.6%) 81 (6.9%)	1.27 (0.95-1.70)	0.54
	yes	18 (4.9%) 19 (4.8%)	1.02 (0.53-1.96)	
Overall		122 (7.8%) 100 (6.3%)	1.23 (0.94-1.61)	

.25 .5 1 2 4

Favors BioFreedom Favors Orsiro

Principal investigator

Prof. Lisette Okkels Jensen, Odense University Hospital, Odense, Denmark

1. Okkels L et al. A Randomized Trial Comparing a Polymer-Free Coronary Drug-Eluting Stent With an Ultra-Thin Strut Bioresorbable Polymer-Based Drug-Eluting Stent in an All-Coroner Patient Population; Presentation; Presented at: TCT 2018; September, 2018; San Diego, USA; Corrected slides, published online on tctMD, Nov 5, 2018; ClinicalTrials.gov: NCT02623140; 2. Okkels L et al. 2-year outcomes of the randomized SORT OUT IX trial, Polymer-free biolimus stent versus the ultrathin strut biodegradable polymer sirolimus-eluting stent in an all-comers population treated with percutaneous coronary intervention, Presented at euroPCR 2021, May 2021, Paris, FRANCE, ClinicalTrials.gov NCT02623140.

BioFreedom is a trademark or registered trademark of the Biosensors International Group. Orsiro and Orsiro Mission are trademarks or registered trademarks of the BIOTRONIK Group of Companies. Clinical data conducted with Orsiro, Orsiro Mission's predecessor device can be used to illustrate Orsiro Mission clinical outcomes.

Manufacturer:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(800) 291-0470 (fax)
www.biotronik.com

Rx ONLY
©2021 BIOTRONIK, Inc.
All rights reserved.
MI419 11/30/21

BIOTRONIK
excellence for life

