




The next level of deliverability¹



Ultrathin struts^{2,a}



Outstanding patient outcomes³




Technical data / ordering info

Vascular Intervention // **Coronary**
Drug-Eluting Stent System



Orsiro[®] Mission DES

Even better deliverability for the outstanding Orsiro DES




Orsiro Mission^{DES}

The next level of deliverability¹


1st in Push⁴

Transmits more force from hub to tip.






1st in Track⁴

Less force needed to follow the path.



1st in Cross⁴


Less force needed to successfully cross demanding anatomies.





proBIO™


Passive coating
for high
biocompatibility



NEW
Deep embedding
for high cross

BIOlute®


Bioabsorbable coating
with controlled
drug release and low
thrombogenicity⁵




Ultrathin
60 μm^a struts
for early
endothelialization



NEW
More flexible shaft
for high track



Enhanced force transmission
for high push



Dual-coating
on shaft for
limited friction



NEW
Ergonomic hub
with kink resistance



^a Nominal strut thickness for size
ø 2.25 - 3.0 mm, mean diameter 62 μm.






Orsiro Mission^{DES}

Ultrathin struts^{2,a} –
thinnest available in
the US⁶

Thinner struts make the difference⁷


- Less disrupted flow
- Improved re-endothelialization

Strut thickness in perspective⁸




For early endothelialization

Strut coverage¹¹
30 days^b




>80%
589 struts analyzed

Strut coverage¹¹
90 days^b



>97%
874 struts analyzed

Strut coverage¹¹
180 days^b




>98%
1,130 struts analyzed



BIO-RESORT Small Vessels (n = 1,506)

Low stent thrombosis (ST) at 5 years¹²



50%
lower ST rate
with **Orsiro** DES^c
p = 0.22

^a Nominal strut thickness for size \varnothing 2.25 - 3.0 mm, mean diameter 62 μm .

^b Images: Secco G et al. Time-related changes in neointimal tissue coverage following a new generation. SES implantation: an OCT observational study. Presented at: euro PCR, May 20, 2014; Paris, France.


^c In comparison to Resolute Integrity, based on 5-year results of the BIO-Resort trial SV subgroup.



Orsiro Mission^{DES}

Outstanding patient outcomes³

One of the most studied DES^d




STUDY NAME	STUDY TYPE	PATIENTS	STATUS	PRIMARY ENDPOINT
BIOSTEMI	RCT	1,300	24-month FU available	TLF at 12 months
TAGLIERI et al.	Network Meta-Analysis	99,039	-	TLF at 12 months and the longest FU available
BIOFLOW-V	RCT	1,334	Completed 60-month FU available	TLF at 12 months
BIO-RESORT	RCT	3,514	Completed 60-month FU available	TVF at 12 months
BIONYX	RCT	2,488	36-month FU available	TVF at 12 months
BIO-SCIENCE	RCT	2,119	Completed, 60-month FU available	TLF at 12 months

Taglieri et al. network meta-analysis

Orsiro DES – the highest probability (70.8%) to rank as the best stent^{14,e}

99,039 patients in a network meta-analysis of 77 RCTs¹⁴

Orsiro DES is associated with a lower 1-year rate of TLF compared with Xience (OR [95% CI] 0.84 [0.71, 0.98], p = 0.03) and Resolute^f (OR [95% CI] 0.81 [0.68, 0.95], p = 0.01).





BIOFLOW-V trial



Pushing the boundaries of performance with Orsiro DES¹⁵

BIOFLOW-V (n = 1,334) FDA pivotal trial

TLF and components out to 5 years^{16,17,18}



ST events out to 5 years¹⁶



TLF – Target Lesion Failure; TV-MI – Target Vessel Myocardial Infarction; TLR – Target Lesion Revascularization; ST – Stent Thrombosis.

^d In large RCTs, based on Taglieri et al. Meta-analysis, against currently used DES.

^e Based on 1-year TLF SUCRA score, in comparison to Xience, Resolute and Nobori/BioMatrix, after a median follow-up period of 50 months.¹⁴

^f Resolute Integrity and Resolute Onyx.

^g The Nobori, BioMatrix, Cre8, Biofreedom and Cypher drug eluting stents are not available in the US.

^h p-values for 60-m frequentist analysis of BIOFLOW-V.¹⁶


ⁱ In comparison to Xience, based on 60-m frequentist analysis of BIOFLOW-V.¹⁶



Orsiro[®] Mission DES

Sirolimus-Eluting Coronary Stent System

Vascular
Intervention
Coronary



Indication

Orsiro Mission is indicated for improving coronary luminal diameter in patients, including those with diabetes mellitus, with symptomatic heart disease, stable angina, unstable angina, non-ST-elevation myocardial infarction or documented silent ischemia due to atherosclerotic lesions in the native coronary arteries with a reference vessel diameter of 2.25 mm to 4.0 mm and a lesion length of \leq 36 mm.

Technical Data

Stent

Stent material	Cobalt chromium, L-605
Strut thickness	\varnothing 2.25 – 3.0 mm: 60 μm^{a} (0.0024"); \varnothing 3.50 – 4.0 mm: 80 μm (0.0031")
Passive coating	proBIO[™] (Amorphous Silicon Carbide)
Active coating	BIOlute[®] bioabsorbable drug matrix consisting of sirolimus and polymer poly-L-lactide (PLLA)
Drug dose	1.4 $\mu\text{g}/\text{mm}^2$

Delivery system

Catheter type	Fast-exchange
Recommended guide catheter	5F (min. I.D. ^j 0.056")
Guide wire diameter	0.014" (0.36 mm)
Usable catheter length	140 cm
Balloon material	Semi crystalline polymer
Coating (distal shaft)	Hydrophilic
Coating (proximal shaft)	Hydrophobic
Marker bands	Two platinum-iridium markers
Proximal shaft diameter	2.0 F
Distal shaft diameter	2.7 F: \varnothing 2.25 – 3.0 mm; 2.9 F: \varnothing 3.5 – 4.0 mm
Nominal pressure (NP)	10 atm
Rated burst pressure (RBP)	16 atm

Ordering Information

Stent \varnothing (mm)	Stent Length (mm)								
	9	13	15	18	22	26	30	35	40
2.25	453925	453931	453937	453943	453949	453955	453961		
2.5	453926	453932	453938	453944	453950	453956	453962	453968	453974
2.75	453927	453933	453939	453945	453951	453957	453963	453969	453975
3.0	453928	453934	453940	453946	453952	453958	453964	453970	453976
3.5	453929	453935	453941	453947	453953	453959	453965	453971	453977
4.0	453930	453936	453942	453948	453954	453960	453966	453972	453978

^a Nominal strut thickness for size \varnothing 2.25 - 3.0 mm, mean diameter 62 μm ; ^j I.D. = Inner Diameter. 1. In comparison to Xience Sierra, Resolute Onyx and Synergy for bench tests on pushability, trackability and crossability, BIOTRONIK data on file; 2. As characterized with respect to strut thickness in Bangalore et al. Newer-Generation Ultrathin Strut Drug-Eluting Stents Versus Older Second-Generation Thicker Strut Drug-Eluting Stents for Coronary Artery Disease: Meta-Analysis of Randomized Trials. *Circulation*. 2018 Nov 13;138(20):2216-26.; 3. Based on investigator's interpretation of BIOFLOW-V primary endpoint result. *Lancet*. 2017 Oct 21; 390(10105):1843-1852; 4. BIOTRONIK data on file; 5. Per investigators' interpretation of preclinical studies with Orsiro as mentioned in Cassese et al. *J Thorac Dis* 2018;10(2):688-692; 6. When compared to FDA approved Drug Eluting Stents. BIOTRONIK data on file; 7. Foin N et al. *International journal of cardiology*. 2014 Dec 20;177(3):800-8; 8. Stefanini GG et al. *Coronary stents: novel developments*. *Heart*. 2014 Jul 1;100(13):1051-61; 9. Low AF. *Stent platform for procedural success: Introducing the Continuous Sinusoidal & Core Wire Technologies*. Presented at: AsiaPCR; 22-24 January, 2015; Singapore, Singapore; 10. Tolentino A. *Evolving DES Strategy: Biodegradable Polymer vs. Bioabsorbable Scaffold*. Presented at: Cardiovascular Nurse/Technologist Symposium; June 17, 2016; New York, USA; 11. Secco G et al. *Time-related changes in neointimal tissue coverage of a novel Sirolimus eluting stent: Serial observations with optical coherence tomography*. *Cardiovascular Revascularization Medicine* 2016; 17(1): 38-43; 12. Ploumen et al. *BIO-RESORT Small Vessels 5Y-EuroPCR2022*; 13. BIOTRONIK data on file, as of January 2020; 14. Taglieri N et al. *Target lesion failure with current drug-eluting stents: Evidence from a comprehensive network meta-analysis*. *JACC* 2020 13(24):2868-78; 15. In comparison to Xience, based on statistically significant lower TV-MI and late/very late definite/probable ST rates from the BIOFLOW-V trial through 5 years; 16. Kandzari D *Ultrathin Bioresorbable Polymer Sirolimus-Eluting Stents versus Thin Durable Polymer Everolimus-Eluting Stents for Coronary Revascularization: Final 5-year Outcomes from the Randomized BIOFLOW V Trial*. Presented at CRT 2022; 17. Kandzari D, et al. *BIOFLOW-V: A Prospective Randomized Multicenter Study to Assess the Safety and Effectiveness of the Orsiro Sirolimus Eluting Coronary Stent System in the Treatment Of Subjects With up to Three De Novo or Restenotic Coronary Artery Lesions*. Science. Presentation at ESC 2017; 18. Kandzari D et al. *Ultrathin, bioresorbable polymer sirolimus-eluting stents versus thin, durable polymer everolimus-eluting stents in patients undergoing coronary revascularisation (BIOFLOW V): a randomised trial*. *Lancet*. 2017 Oct 21; 390(10105):1843-1852.

Orsiro, Orsiro Mission, **proBIO** and **BIOlute** are trademarks or registered trademarks of the BIOTRONIK Group of Companies. Synergy, Synergy XD, Promus and Taxus are trademarks or registered trademarks of the Boston Scientific Group of Companies. Resolute, Integrity, Resolute Onyx and Resolute Integrity are trademarks or registered trademarks of the Medtronic Group of Companies. Xience and Xience Skypoint are trademarks or registered trademarks of the Abbott Group of Companies. Biofreedom and BioMatrix are trademarks or registered trademarks of the Biosensors International Group. Cre8 is a trademark or registered trademark of the Alvimedica Group of Companies. Nobori is a trademark or registered trademark of the Terumo Group of Companies. Cypher is a trademark or registered trademark of the Cardinal Health Group of Companies.

Manufacturer:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com



Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
Tel (800) 547-0394 (24-hour)
Fax (800) 291-0470
www.biotronik.com

©2022 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification,
revision and improvement.

BIOTRONIK
excellence for life

